

SMALL DUAL-SENSOR LOGGERS



COMPACT,
ACCURATE,
DEPENDABLE

RBRduet⁴|2x loggers are compact, lightweight instruments designed for high-precision oceanographic measurements in demanding environments. They combine exceptional accuracy and resolution using any standard AA battery chemistry. Engineered for durability in harsh marine conditions, the RBRduet⁴|2x delivers reliable, long-term performance across a wide range of oceanographic applications.


FEATURES



Compact and
lightweight



Two AA
batteries



Up to 32Hz
sampling



USB-C
download



Infinite
memory*

*not really, but we stopped counting at billions of samples.

Available configurations

- ▶ RBRduet⁴ D.ODO|2x pressure and dissolved oxygen; depths up to 1000m
- ▶ RBRduet⁴ D.Tu|2x pressure and turbidity; depths up to 1000m
- ▶ RBRduet⁴ D.chl-a|2x pressure and chlorophyll-a; depths up to 1000m

Deep variant

- ▶ |deep depths up to 6000m

SMALL DUAL-SENSOR LOGGERS

COMPACT, ACCURATE, DEPENDABLE

Specifications

Physical

Configuration	D.Tu	D.chl-a	D.ODO
Storage	Infinite memory ¹		
Power	Any AA cells, any chemistry		
Communication	USB-C		
Clock drift	±60 seconds per year		
Housing	Plastic (pl) or Titanium (Ti)		
Max depth rating	1700m (pl) ² 10000m (Ti)		
Diameter	~33mm		
Length	374mm (pl) 379mm (Ti)	425mm (pl) 425mm (Ti)	
Weight in air in water	254g (pl) 560g (Ti) 55g (pl) 360g (Ti)	304g (pl) 625g (Ti) 65g (pl) 390g (Ti)	

¹Not really, but we stopped counting at billions of samples.

²Actual depth rating is determined by pressure sensor.

Dissolved oxygen

Measurement range	0-1000µmol/L
Calibrated range	0-500µmol/L concentration 0 – 120% saturation 1.5°C to 30°C temperature
Initial accuracy	Max of ±8µmol/L or ±5% fast Max of ±2µmol/L or ±1.5% standard Max of ±2µmol/L or ±1.5% slow
Resolution	<1µmol/L (saturation 0.4%) fast <0.5µmol/L (saturation 0.2%) standard <0.1µmol/L (saturation 0.04%) slow
Time constant	<1s fast, <8s standard, or <30s slow
Sampling rates	24hr to 1Hz
Output Values	Temperature (°C) Dissolved O ₂ concentration (µmol/L) Dissolved O ₂ concentration (salinity compensated, µmol/L) Dissolved O ₂ saturation (%) Dissolved O ₂ phase (°)

Pressure

Range	20 / 50 / 100 / 200 / 500 / 1000dbar (plastic) 1000 / 2000 / 6000dbar (Ti)
Initial accuracy	±0.05% full scale
Resolution	<0.001% full scale
Typical stability	±0.05% full scale / year
Time constant	<10ms

RBR Ltd

+1 613 599 8900
info@rbr-global.com
rbr-global.com

Temperature

Range	-5°C to 35°C
Initial accuracy	±0.002°C
Resolution	<0.00005°C
Typical stability	±0.002°C / year
Time constant	<1s

Chlorophyll-a

Wavelength	470nm/695nm (excitation/emission)
Calibrated range*	0-50µg/L
Measurement range*	0-500µg/L
Detection limit	0.020µg/L
Optical	
Linearity, R ²	0.99
Initial accuracy	5%

*Scaled for the in vivo fluorescence response

Turbidity

Wavelength	880nm
Centroid angle	90°
Linearity, R ²	0.99
Initial accuracy	5%
Calibrated range	0 - 1000FTU
Measurement range	0 - 1500FTU ¹
Detection limit	0.005FTU
Optical backscatter	
Wavelength	880nm
Centroid angle	135°
Linearity, R ²	0.99
Initial accuracy	5%
Calibrated range	1000 - 4000FTU
Measurement range	0 - 20000FTU ²
Detection limit	2.0FTU

¹ Response becomes non-linear above 1000FTU.

² Response becomes non-linear below 500FTU and above 15000FTU.

